


# Float Assembly - Model 30 & 46

1. Float Bulb - 200082
2. Float Stem - 400001
3. Teflon Wrapped Float Guide - 400008
4. Float Gasket - Silicone -----
5. Float Diaphragm-----400277
6. Washer Stainless Steel - 300177
7. Float Switch Actuator - 200052
8. Stainless Steel Float Washer - 100024
9. Stainless Steel 1/4 x 20 Nut - 300125
10. Stainless Steel Washers (1/4 FCW) - 300176



# DURASTILL

4200 NE Birmingham Rd. • Kansas City, MO 64117  
(816) 454-5260 • FAX (816) 452-7581